## Photometer AL100

Molybdän • Molybdenum •  
Molybdène • Molibdeno • Molibdato

Tablet • Powder

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Seite 3–19

**(GB)** Instruction Manual

Page 21–37

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Page 39–55

**(IT)** Istruzioni d'uso

Pagina 57–73

**(ES)** Instrucciones

Página 75–91

**CE-Konformitätserklärung / Declaration of CE-Conformity  
Déclaration de conformité CE / Dichiarazione di conformità CE /  
CE-Declaración de conformidad**

---

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
• <b>General notes</b> .....	22
Guidelines for photometric measurements .....	22
Method notes .....	22
Replacement of batteries .....	23
• <b>Functional description</b> .....	24
Operation .....	24
Display backlight .....	25
Recall of stored data .....	25
Countdown .....	25
• <b>Methods</b> .....	26
Molybdenum LR, with powder pack (Mo 1) .....	26
Molybdenum HR, with powder pack (Mo 2) .....	28
Molybdenum, with tablet (Mo 3) .....	30
• <b>Menu options</b> .....	32
Menu selections .....	32
Recall of stored data .....	32
Transmitting stored data .....	32
Setting date and time .....	33
• <b>Calibration Mode</b> .....	33
User calibration .....	33
Factory calibration reset .....	35
• <b>Technical data</b> .....	36
Operating messages .....	37
Error codes .....	37

**CAUTION**

The accuracy of the instrument is only valid if the instrument is used in an environment with controlled electromagnetic disturbances according to DIN 61326. Wireless devices, e.g. wireless phones, must not be used near the instrument.


## Guidelines for photometric measurements

1. Vials, caps and stirring rods should be cleaned thoroughly **after each analysis** to prevent interference. Even minor reagent residues can cause errors in the test result.
2. The outside of the vial must be clean and dry before starting the analysis. Clean the outside of the vials with a towel to remove fingerprints or other marks.
3. Zero calibration and test must be carried out with the same vial as there may be slight differences in optical performance between vials.
4. The vials must be positioned in the sample chamber for zeroing and test with the  $\Delta$  mark on the vial aligned with the  $\nabla$  mark on the instrument.
5. Always perform zeroing and test with the vial cap tightly closed. Only use the cap with a sealing ring.
6. Bubbles on the inside wall of the vial lead to incorrect measurements. To prevent this, remove the bubbles by swirling the vial before performing the test.
7. Avoid spillage of water into the sample chamber because this can lead to incorrect test results.
8. Contamination of the transparent cell chamber can result in wrong readings. Check at regular intervals and – if necessary – clean the transparent cell chamber using a moist cloth or cotton buds.
9. Large temperature differences between the instrument and the environment can lead to errors – e.g. due to the formation of condensation in the cell chamber or on the vial.
10. To avoid errors caused by stray light do not use the instrument in bright sunlight.
11. Always add the reagent tablets to the water sample straight from the foil without touching them with the fingers.
12. The reagents must be added in the correct sequence.


## Method notes

- Prior to measurement ensure that the sample is suitable for analysis (no major interferences) and does not require any preparation i.e. pH adjustment, filtration etc.
- Reagents are designed for use in chemical analysis only and should be kept well out of the reach of children.
- Ensure proper disposal of reagent solutions.
- Material Safety Data Sheets are available on request (Internet: [www.aqualytic.de](http://www.aqualytic.de))


**Correct position of the vial (Ø 24 mm):**

**Correct filling of the vial:**

**Replacement of batteries:**

**CAUTION:**

**To ensure that the instrument is water proof:**

- seal ring (E) must be in position
- battery compartment cover (B) must be fixed with the four screws

If the batteries are removed for more than one minute the date and time menu starts automatically when the photometer is switched on the next time.

## Operation

METHOD

Switch the unit on using the [ON/OFF] key.

The display shows the following:

Select the required test using the [MODE] key.

### Scroll Memory (SM)

To avoid unnecessary scrolling for the required test method, the instrument memorizes the last method used before being switched off. When the instrument is switched on again, the scroll list comes up with the last used test method first.

METHOD

The display shows the following:

Fill a clean vial with the water sample up to the 10 ml mark, screw the cap on and place the vial in the sample chamber making sure that the  $\Sigma$  marks are aligned.

METHOD

0.0.0

Press the [ZERO/TEST] key.

The "Method" symbol flashes for approx. 8 seconds.

The display shows the following:

After zero calibration is completed, remove the vial from the sample chamber. The characteristic coloration appears after the addition of the reagents.

Replace the cap on the vial and place in the sample chamber making sure that the  $\Sigma$  marks are aligned.

METHOD

RESULT

Press the [ZERO/TEST] key.

(For Countdown/reaction period see page 25)

The "Method" symbol flashes for approx. 3 seconds.

The result appears in the display.

The result is saved automatically.

### Repeating the test:

Press the [ZERO/TEST] key again.

### Repeating the zero:

Press the [ZERO/TEST] key for 2 seconds.

### Display backlight

Press the [!] key to turn the display backlight on or off. The backlight is switched off automatically during the measurement.

### Recall of stored data

If the instrument is switched on, press the [!] key for more than 4 seconds to access the recall menu.

### Countdown / reaction period

If a reaction period is included in a method a countdown function can be used:

Press the [!] key and hold.

Press the [ZERO/TEST] key.

Release the [!] key; the countdown starts.

After the countdown is finished the measurement starts automatically.

It is possible to interrupt the countdown by pressing the [ZERO/TEST] key. Measurement starts immediately.

**Caution:**

**An incomplete reaction period can lead to incorrect test results.**

## Molybdenum LR with Vario Powder Pack

### 0.03 – 3.0 mg/l Mo

Fill a clean Mixing Cylinder (25-ml) with **20 ml of the water sample**.

Add the contents of **one Vario Molybdenum 1 LR F20 Powder Pack** straight from the foil into the water sample (20 ml).

Close the Mixing Cylinder tightly with a stopper and swirl several times to dissolve the powder.

Use two clean vials (24 mm Ø) and mark one as blank for zeroing.

Fill each vial with 10 ml of pre prepared water sample.

Close the blank tightly with the cap.

Add **0,5 ml of Vario Molybdenum 2 LR solution** to the sample.

Close the vial tightly with the cap and invert several times to mix the contents.

**Wait for a reaction period of 2 minutes.**

Place the blank in the sample chamber making sure that the  $\Sigma$  marks are aligned.

Press [ZERO/TEST] key.

The method symbol flashes for approx. 8 seconds.



Remove the vial from the sample chamber.

Place the sample in the sample chamber making sure that the  $\Sigma$  marks are aligned.

Press [ZERO/TEST] key.

The method symbol flashes for approx. 3 seconds.

The result is shown in the display in mg/l Molybdenum.

#### Tolerance:

0 – 1 mg/l:  $\pm 0.05$  mg/l Mo

> 1 – 3 mg/l:  $\pm 0.10$  mg/l Mo

#### Conversion:

mg/l  $\text{MoO}_4 = \text{mg/l Mo} \times 1.67$

mg/l  $\text{Na}_2\text{MoO}_6 = \text{mg/l Mo} \times 2.15$

**Notes:**

1. Strong alkaline or acidic water samples must be adjusted between pH 3 and pH 5 before the reagent is added (use 0.5 mol/l Sulfuric acid resp. 1 mol/l Sodium hydroxide).
2. Before using clean the vials with Hydrochloric acid (approx. 20%). Rinse thoroughly with deionised water.

Mo 2

Tablet

Powder

## Molybdenum HR with Vario Powder Pack 0.3 – 40 mg/l Mo

0.0.0

Tablet

Powder

Fill a clean vial (24 mm Ø) with **10 ml of the water sample** and perform zero calibration (see "Operation").

Add the contents of **one Molybdenum HR 1 F10 Vario Powder Pack** straight from the foil into the water sample.

Close the vial tightly with the cap and swirl several times to dissolve the powder.

Add the contents of **one Vario Molybdenum HR 2 F10 Powder Pack** straight from the foil to the same water sample.

Close the vial tightly with the cap and invert several times to mix the contents.

Add the contents of **one Vario Molybdenum HR 3 F10 Powder Pack** straight from the foil to the same water sample.

Close the vial tightly with the cap and swirl several times to dissolve the powder.

Place the vial in the sample chamber making sure that the  $\Sigma$  marks are aligned.

**Wait for a reaction period of 5 minutes.**

(Countdown can be activated, see page 25)

The method symbol flashes for approx. 3 seconds.

The result is shown in the display in mg/l Molybdenum.

Mo 2

Tablet

Powder

RESULT

Tablet

Powder

### Tolerance:

0 – 30 mg/l:  $\pm 2$  mg/l Mo

> 30 – 40 mg/l:  $\pm 3$  mg/l Mo

### Conversion:

mg/l  $\text{MoO}_4$  = mg/l Mo x 1.67

mg/l  $\text{Na}_2\text{MoO}_6$  = mg/l Mo x 2.15

**Notes:**

1. Filter turbid water samples using filter paper and funnel before analysis.
2. Highly buffered water samples or extreme pH values should be adjusted to a pH of nearly 7 with 1 mol/l Nitric acid or 1 mol/l Sodium hydroxide.
3. Concentrations above 10 mg/l Cu causes too high test values if the reaction time of 5 minutes is increased. So it is very important to perform the test procedure as described.
4. Substances which may interfere when present in concentrations at:

Aluminium	50 mg/l
Chromium	1000 mg/l
Iron	50 mg/l
Nickel	50 mg/l
Nitrite	all levels

## Molybdate with Tablet

### 0.6 – 30 mg/l Mo

Fill a clean vial (24 mm Ø) with **10 ml of the water sample** and perform zero calibration (see "Operation").

Remove the vial from the sample chamber and **empty the vial**.

Fill **20 ml of the water sample** in a 100 ml beaker.

Add **one MOLYBDATE HR No. 1 tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.

Add **one MOLYBDATE HR No. 2 tablet** straight from the foil to the same water sample and crush the tablet using a clean stirring rod.

Dissolve the tablets using a clean stirring rod.

Rinse out the vial with the prepared water sample and then fill to the 10 ml mark.


Close the vial tightly with the cap.

Place the vial in the sample chamber making sure that the  $\Sigma$  marks are aligned.

Press the [ZERO/TEST] key.

The method symbol flashes for approx. 3 seconds.

The result is shown in the display in mg/l Molybdenum.

#### Tolerance:

0 – 15 mg/l:  $\pm 1.0$  mg/l Mo  
> 15 – 30 mg/l:  $\pm 1.5$  mg/l Mo

#### Conversion:

$\text{mg/l MoO}_4 = \text{mg/l Mo} \times 1.67$   
 $\text{mg/l Na}_2\text{MoO}_6 = \text{mg/l Mo} \times 2.15$

**Notes:**

1. The tablets must be added in the correct sequence.
2. Under test conditions (pH 3.8 – 3.9) iron does not interfere nor do other metals at levels likely to be found in industrial water systems.

## Menu selections

Press the [MODE] key and **hold**.

Switch the unit on using the [ON/OFF] key.

Allow the 3 decimal points to be displayed before releasing the [MODE] key.

The [!] key allows for selection of the following menu points:

- ▲ diS recall stored data
- ▲ Prt printing stored data
- ▲ setting the date and time
- Cal user calibration

### ▲ diS – Recall of stored data

After confirming the selection with the [MODE] key the photometer shows the last 16 data sets in the following format (automatically proceeds every 3 seconds until result is displayed):

```
Number  n xx (xx: 16...1)
Year    YYYY (e.g. 2010)
Date    mm.dd (monthmonth:dayday)
Time    hh:mm (hourhour:minute)
Test    Method
Result  x,xx
```

The [ZERO/TEST] key repeats the current data set.


The [MODE] key scrolls through all stored data sets.

Quit the menu by pressing [!] key.

### ▲ Prt – Transmitting stored data (to Printer or PC)

**Note:** To print data, or to transmit to a PC, the optional IRiM (Infrared Interface Module) is required.

The IRiM Module and the connected printer/PC must be ready. Press the [MODE] key to start the transmitting, the instrument displays "PrtG" (Printing) for approx. 1 second followed by the number of the first data set and its transmission. All data sets will be transmitted one after the other. After finishing the instrument switches to test mode.

The print job can be cancelled by pressing the [On/Off] key. The instrument switches off.

**E 132**

If the instrument is not able to communicate with the IRiM, a timeout occurs after approx. 2 minutes. The error E 132 is displayed for approx. 4 seconds. Subsequently, the instrument switches to test mode (see also IRiM manual).

**SET**

**DATE**

**YYYY**

(2 sec.)

**⚠ Setting date and time (24-hour-format)**

After confirming the selection with the [MODE] key the value to be edited will be shown for 2 sec.

The setting starts with the year (YYYY) followed by the actual value to be edited. The same applies for month (mm), day (dd), hour (hh) and minutes (mm). Set the minutes first in steps of 10, press the [!] key to continue setting the minutes in steps of 1.

Increase the value by pressing the [MODE] key.

Decrease the value by pressing [ZERO/TEST] key.

Proceed to the next value to be edited by pressing [!] key.

After setting the minutes and pressing the [!] key the display will show "IS SET" and the instrument returns to the measurement mode.

**CAL**

**CAL**

**CAL**

**METHOD**

≡METHOD≡

**0.0.0**

**CAL**

≡METHOD≡

**Cal User calibration**

**Note:**

user calibration (Display in calibration mode)

factory calibration (Display in calibration mode)

After confirming the selection with the [MODE] key the instrument will show CAL"Method".

Scroll through methods using the [MODE] key.

Fill a clean vial with the standard up to the 10 ml mark, screw the cap on and place the vial in the sample chamber making sure that the X marks are aligned.

Press the [ZERO/TEST] key.

The method symbol flashes for approx. 8 seconds.

The display shows the following in alternating mode:

Perform calibration with a standard of known concentration (see "Operation").

Press the [ZERO/TEST] key.

The method symbol flashes for approx. 3 seconds.

## **GB** Calibration Mode

RESULT

CAL

The result is shown in the display, alternating with CAL.

If the reading corresponds with the value of the calibration standard (within the specified tolerance), exit calibration mode by pressing the [ON/OFF] key.

Changing the displayed value:

Mode

Pressing the [MODE] key once increases the displayed value by 1 digit.

Zero  
Test

Pressing the [ZERO/TEST] key once decreases the displayed value by 1 digit.

CAL

Press the corresponding key until the reading equals the value of the calibration standard.

RESULT + x

On  
Off

By pressing the [ON/OFF] key, the new correction factor is calculated and stored in the user calibration software.

: :

Confirmation of calibration (3 seconds).

## Factory calibration reset

Resetting the user calibration to the original factory calibration will reset all methods and ranges.

A user calibrated method is indicated by a "Cal" symbol while the test result is displayed.

To reset the calibration press both the [MODE] and [ZERO/TEST] key and **hold**.

Switch the unit on using the [ON/OFF] key.

Release the [MODE] and [ZERO/TEST] keys after approx. 1 second.

The following messages will appear in turn on the display:


The factory setting is active.  
(SEL stands for Select)

**or:**

Calibration has been set by the user.

(If the user calibration is to be retained, switch the unit off using the [ON/OFF] key).

Calibration is reset to the factory setting by pressing the [MODE] key.

The following messages will appear in turn on the display:

Switch the unit off using the [ON/OFF] key.

**Technical Data**

Instrument	double wavelength, automatic wavelength selection, direct reading colorimeter
Light source:	LEDs, interference filters (IF) and photosensor in transparent cell chamber. Wavelength specifications of the IF: 430 nm $\Delta \lambda = 5 \text{ nm}$ 610 nm $\Delta \lambda = 6 \text{ nm}$
Wavelength accuracy	$\pm 1 \text{ nm}$
Photometric accuracy*	3% FS (T = 20° C – 25° C)
Photometric resolution	0.01 A
Power supply	4 batteries (AAA/LR 03) lifetime: approx. 17 hours or 5000 tests
Auto-OFF	automatic switch off 15 minutes after last keypress
Display	backlit LCD (on keypress)
Storage	internal ring memory for 16 data sets
Interface	IR interface for data transfer
Time	real time clock und date
Calibration	user and factory calibration resetting to factory calibration possible
Dimensions	155 x 75 x 35 mm (LxWxH)
Weight	approx. 260 g (incl. batteries)
Ambient conditions	temperature: 5–40 °C rel. humidity: 30–90% (non-condensing)
Waterproof	floating; as defined in IP 68 (1 hour at 0.1 meter)
CE	Certificate for Declaration of CE-Conformity at <a href="http://www.aqualytic.de">www.aqualytic.de</a>

*\*measured with standard solutions*

To ensure maximum accuracy of test results, always use the reagent systems supplied by the instrument manufacturer.

## Operating messages

**Hi**

Measuring range exceeded or excessive turbidity.

**Lo**

Result below the lowest limit of the measuring range.

Replace batteries, no further tests possible.

**btLo**

Battery capacity is too low for the display backlight; measurement is still possible.

Store  
Print  
Cal  
Date  
Time  
**RESULT**

A user calibrated method is indicated by a "Cal" symbol while the test result is displayed. (see "Factory calibration reset").

## Error codes

**E27 / E28 / E29**

Light absorption too great. Reasons: e.g. dirty optics.

**E 10 / E 11**

Calibration factor "out of range"

**E 20 / E 21**

Too much light reaching the detector.

**E23 / E24 / E25**

Too much light reaching the detector.

**E 22**

Battery capacity was too low during measurement. Change battery.

**E 70**

Mo 1 (Powder Pack, LR): Factory calibration incorrect / erased

**E 71**

Mo 1 (Powder Pack, LR): User calibration incorrect / erased

**E 72**

Mo 2 (Powder Pack, HR): Factory calibration incorrect / erased

**E 73**


Mo 2 (Powder Pack, HR): User calibration incorrect / erased

**E 74**

Mo 3 (Tablet): Factory calibration incorrect / erased

**E 75**

Mo 3 (Tablet): User calibration incorrect / erased

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**Tintometer GmbH, Division Aqualytic®**